


# NEW HEARTMATE 3™

*with Full MagLev™ flow technology*

RESPECTS  
THE FLOW



 THORATEC®

HeartMate 3™

with Full MagLev™

*New HeartMate 3™ with Full MagLev™ flow technology*

## HAEMOCOMPATIBILITY MAKES THE DIFFERENCE



Full MagLev flow technology is designed to optimise haemocompatibility and reduce blood trauma through gentle blood handling.



## Remarkably low incidence of key adverse events

**0**  
Pump thrombosis<sup>1</sup>

**0**  
Haemolysis<sup>1</sup>

**0**  
Pump malfunctions or exchanges<sup>1</sup>

## Excellent 30-day and 6-month survival

**98%**  
30-day survival<sup>1</sup>

**92%**  
6-month survival<sup>1</sup>

**Study description:** Single-arm, prospective, multicentre, nonblinded, nonrandomised study to evaluate the performance and safety of HeartMate 3 Left Ventricular Assist System after 6 months of support. Patients were adults with an ejection fraction  $\leq 25\%$  and cardiac index  $\leq 2.2$  L/min/m<sup>2</sup> while not on inotropes. Patients must have also met one of the following: on optimal medical management, based on current heart failure practice guidelines, for at least 45 out of the previous 60 days and were failing to respond; or in NYHA class III or IV heart failure for at least 14 days AND dependent on intra-aortic balloon pump for at least 7 days; or inotrope-dependent/unable to wean from inotropes; or listed for transplant.<sup>1</sup>


**Patient population:** Fifty patients were enrolled and implanted with HeartMate 3 at 10 centres. Patients were representative of an advanced heart failure population and were a mix of bridge-to-transplant (54%) and destination therapy (46%). Twenty-six patients (52%) were in INTERMACS profiles 2 and 3, and 24 patients (48%) were in profiles 4 to 6. No patients were in profiles 1 or 7.

**HeartMate 3™**  
with Full MagLev™

# GENTLE BLOOD HANDLING

with Full MagLev™ flow technology


DESIGNED FOR HAEMOCOMPATIBILITY  
with the goal of minimising complications



- Rotor is completely suspended by magnetic forces, preventing surface-to-surface contact that could cause blood trauma
- Designed for:
  - Low haemolysis
  - Minimised blood shearing and stasis

## Large, consistent blood-flow pathways

Blood-flow gaps 10 to 20 times greater than those of hydrodynamic-bearing pumps



Dimensions are approximate for illustration purposes only.

The use of a red blood cell as a measuring unit is for illustration purposes and is not meant to imply actual blood-flow quantities during operation.

## Full range of operation

Provides flow from 2.5 to 10 L/min to accommodate a broad range of clinical needs


## Artificial pulsatility

Designed to promote pump washing and minimise areas of stasis in the pump

# PERFORMANCE-LEADING OUTCOMES


with HeartMate 3™

Reach significant improvements in NYHA class<sup>1</sup>



Gain significant improvement in 6-minute walk distance<sup>2</sup>


Average 6-minute walk distance in patients implanted with HeartMate 3 (n = 36)



NYHA: New York Heart Association.

See progressive, sustained QoL improvements<sup>1</sup>

EQ-5D Visual Analog Scale



Mean baseline values for paired analyses were 52.4, 52.1, and 51.4 at 1, 3, and 6 months, respectively.

Confirm a favourable safety profile<sup>1</sup>

Adverse events through 6 months for patients implanted with HeartMate 3 (N = 50)

ADVERSE EVENTS	Number of Patients	Percentage of Patients	Number of Events
Device thrombosis	0	0%	0
Haemolysis	0	0%	0
Stroke*			
Ischaemic	4	8%	4
Haemorrhagic	2	4%	2
Bleeding requiring surgery	7	14%	8
GI bleeding	4	8%	6
Device malfunction	0	0%	0
Driveline infection	5	10%	5
Right heart failure	5	10%	5

\*Includes 3 procedure-related events: 1 implant issue (difficulty engaging inflow conduit); 1 following anaphylactic shock from contrast media; 1 following transcatheter aortic valve implantation procedure.

# Discover a better experience

- Built to the highest standards of Swiss engineering
- Designed for intrapericardial placement
- Features a thin, mechanical apical cuff lock for quick and easy pump attachment
- Incorporates a modular driveline for straightforward replacement of external portion
- Offers an advanced controller and extended battery life for simplified patient management



*HeartMate 3™ treats blood with the respect that it deserves,  
for optimised outcomes and better patient lives.*

0

Pump  
thrombosis<sup>1</sup>

0

Haemolysis<sup>1</sup>

0

Pump malfunctions  
or exchanges<sup>1</sup>

98%

30-day survival<sup>1</sup>

92%

6-month survival<sup>1</sup>

References: 1. Netuka I, et al. HeartMate 3 fully magnetically levitated LVAD for the treatment of advanced heart failure: results from the CE Mark Trial. Presented at 19th Annual Meeting of the Heart Failure Society of America; September 26-29, 2015; Washington, DC. 2. Data on file. Thoratec Corporation, 2015.

© 2015 Thoratec Corporation. All rights reserved. EU-HM3-08150343

Thoratec and the Thoratec logo are registered trademarks and HeartMate 3 and Full MagLev are trademarks of Thoratec Corporation.

Not approved for commercial sale in the US.

 THORATEC®

HeartMate 3™

with Full MagLev™