


Dragonfly™ OPTIS™

Imaging Catheter




PRODUCT FEATURES

- OCT imaging catheter for OPTIS™ Integrated System*, ILUMIEN™ OPTIS™ and ILUMIEN™ PCI Optimization™ Systems
- Radiopaque markers at the distal tip, imaging lens and 50 mm proximal to the lens for aligning the OCT pullback to region of interest
- Proprietary sheath design provides continuous calibration throughout pullback
- Dual lumen rapid exchange tip for easier guidewire loading
- Lens marker <2 mm proximal to optical lens
- Proximal shaft marker bands at 90 mm and 100 mm
- 2.7 F low-profile tip
- Hydrophilic coating
- Simple setup
- Compatible with standard 0.014" guidewires

SPECIFICATIONS

Usable length:	135 cm
Outer diameter:	2.7 F (distal)
Wire lumen:	0.014"



*Available only in geographies with appropriate regulatory approvals

ORDERING INFORMATION

Reorder Number	Description
C408646	Dragonfly OPTIS Imaging Catheter, Sterile DOC Cover and 1 Sterile 3 ml Syringe are provided with the catheter
200-700-00	Sterile DOC Cover
C408655	Sterile 3 ml Syringe

St. Jude Medical Inc.
Global Headquarters
One St. Jude Medical Drive
St. Paul, Minnesota 55117
USA
+1 651 756 2000
+1 651 756 3301 Fax

St. Jude Medical
International Division
One Lillehei Plaza
St. Paul, Minnesota 55117
USA
+1 651 756 2000
+1 651 756 2291 Fax

St. Jude Medical S.C., Inc.
Americas Division
6300 Bee Cave Road
Bldg. Two, Suite 100
Austin, TX 78746
USA
+1 512 286 4000
+1 512 732 2418 Fax

SJM Coordination Center BVBA
The Corporate Village
Da Vincilaan 11-Box F1
B-1935 Zaventem, Belgium
+32 2 774 68 11
+32 2 772 83 84 Fax

St. Jude Medical Brasil Ltda.
Rua Itapeva, 538
5º ao 8º andares
01332-000 – São Paulo – SP
Brazil
+55 11 5080 5400
+55 11 5080 5423 Fax

St. Jude Medical (Hong Kong) Ltd.
Suite 1608, 16/F Exchange Tower
33 Wang Chiu Road
Kowloon Bay, Kowloon
Hong Kong SAR
+852 2996 7688
+852 2956 0622 Fax

St. Jude Medical Japan Co., Ltd.
Shiodome City Center 15F
1-5-2 Higashi Shinbashi, Minato-ku
Tokyo 105-7115
Japan
+81 3 6255 6370
+81 3 6255 6371 Fax

St. Jude Medical Australia Pty, Ltd.
17 Orion Road
Lane Cove, NSW 2066
Australia
+61 2 9936 1200
+61 2 9936 1222 Fax

ILUMIEN.com
SJMprofessional.com

**Rx Only**

Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Product referenced is approved for CE Mark.

Dragonfly OPTIS catheter, OPTIS integrated system and ILUMIEN system

Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL and the nine-squares symbol are trademarks and service marks of St. Jude Medical, Inc. and its related companies. © 2014 St. Jude Medical, Inc. All Rights Reserved.

ID-2001610 A EN (08/14)