

Angio-Seal™ Evolution™

Vascular Closure Device

Product Highlights

- Automated collagen compaction system for standardized deployment
- Standardized deployment provides more control which may reduce procedural variables
- Instant compaction is designed to provide control of the arteriotomy, enabling an optimal seal
- Suture, collagen and anchor are fully bioabsorbable within 60–90 days

Ordering Information

Contents: Vascular Closure Device, Insertion Sheath, Arteriotomy Locator and 70 cm Guidewire with “J” Straightener
(10 units per box)

| Reorder Number | French Size | Guidewire Diameter (in) |
|----------------|-------------|-------------------------|
| C610136 | 6 | 0.035 |
| C610137 | 8 | 0.038 |



Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Product referenced is approved for CE Mark.

Angio-Seal, Evolution, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies.

Angio-Seal™ VIP Vascular Closure Device

Product Highlights

- VIP V-Twist Integrated Platform is a unique technology that provides a larger collagen footprint for greater arteriotomy coverage
- Bondek®-Plus suture with a patented Polyglyd™ coating for maximum lubricity
- Suture, collagen and anchor are fully bioabsorbable within 60–90 days

Ordering Information

Contents: Vascular Closure Device, Insertion Sheath, Arteriotomy Locator and 70 cm Guidewire with “J” Straightener
(10 units per box)

| Reorder Number | French Size | Guidewire Diameter (in) |
|----------------|-------------|-------------------------|
| 610132 | 6 | 0.035 |
| 610133 | 8 | 0.038 |



Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Bondek-Plus is a registered trademark of TFX Medical. Polyglyd is a registered trademark of Teleflex-CT Devices Incorporated. Product referenced is approved for CE Mark.

Angio-Seal, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies.

Angio-Seal™ STS Plus Vascular Closure Device

Product Highlights

- Self-tightening suture provides a secure anchor collagen sandwich
- Secure-Cap™ provides tactile and audible confirmation of device positioning during deployment
- Suture, collagen and anchor are fully bioabsorbable within 60–90 days

Ordering Information

Contents: Vascular Closure Device, Insertion Sheath, Arteriotomy Locator and 70 cm Guidewire with “J” Straightener
(10 units per box)

| Reorder Number | French Size | Guidewire Diameter (in) |
|----------------|-------------|-------------------------|
| 610120 | 6 | 0.035 |
| 610122 | 8 | 0.038 |



Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Product referenced is approved for CE Mark.

Angio-Seal, Secure-Cap, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies.