


St. Jude Medical™ Quartet™ Family of LV Leads




# PROVEN QUARTET™ LV LEAD PERFORMANCE

WITH THE MOST QUADRIPOLEAR LEAD OPTIONS TO MATCH A PATIENT'S ANATOMY

## Built on the proven Quartet™ LV lead platform and featuring:

- Four uniquely spaced electrodes to provide flexibility in pacing location
- Small 4.7 F lead body and 4.0 F lead tip
- Allows patients to safely undergo an MRI scan when used in combination with a St Jude Medical™ MRI Ready device<sup>1,2</sup>
- Optim™ lead insulation




The Quartet™ Family of LV leads offers more distal shape options including the Large-S, Small-S and the Double Bend and more total electrode spacing options including 40, 47 and 60 mm to ensure not only successful CRT Implant but true quadripolar flexibility and options to help patient optimization.




**1457Q “Double-Bend” Quartet™ LV Lead**

- May be ideal for medium size and length target vessels
- Total spacing 47 mm (20-10-17)




**1456Q “Small-S” Quartet™ LV Lead**

- May be ideal for shorter target vessels
- Total spacing 40 mm (20-10-10)




**1458QL “Wide-Spaced” Quartet™ LV Lead**

- May be ideal for larger size and longer length target vessels
- Total spacing 60 mm (20-27-13)



**1458Q “Original Large-S” Quartet™ LV Lead**

- Superb deliverability with exceptional stability and performance
- Total spacing 47 mm (20-10-17)



**Brief Summary:** Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Unless otherwise noted, <sup>TM</sup> indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL and the nine-squares symbol are trademarks and service marks of St. Jude Medical, Inc. and its related companies. © 2016 St. Jude Medical, Inc. All Rights Reserved.

**St. Jude Medical Inc.**

Global Headquarters  
One St. Jude Medical Drive  
St. Paul, MN 55117  
USA  
**T** +1 651 756 2000 | **F** +1 651 756 3301

**St. Jude Medical S.C., Inc.**

Americas Division  
6300 Bee Cave Road  
Bldg. Two, Suite 100  
Austin, TX 78746  
USA  
**T** +1 512 286 4000 | **F** +1 512 732 2418

**SJM Coordination Center BVBA**

The Corporate Village  
Da Vincilaan 11-Box F1  
B-1935 Zaventem, Belgium  
**T** +32 2 774 68 11 | **F** +32 2 772 83 84

**St. Jude Medical Brasil Ltda.**

Rua Itapeva, 538  
5º ao 8º andares  
01332-000 – São Paulo – SP  
Brazil  
**T** +55 11 5080 5400 | **F** +55 11 5080 5423

**St. Jude Medical (Hong Kong) Ltd.**

Suite 1608, 16/F Exchange Tower  
33 Wang Chiu Road  
Kowloon Bay, Kowloon  
Hong Kong SAR  
**T** +852 2996 7688 | **F** +852 2956 0622

**St. Jude Medical Japan Co., Ltd.**

Shiodome City Center 15F  
1-5-2 Higashi Shinbashi, Minato-ku  
Tokyo 105-7115  
Japan  
**T** +81 3 6255 6370 | **F** +81 3 6255 6371

**St. Jude Medical Australia Pty, Ltd.**

17 Orion Road  
Lane Cove, NSW 2066  
Australia  
**T** +61 2 9936 1200 | **F** +61 2 9936 1222

**SJM.com**