

Product Highlights

- Options for anatomic variabilities
- Multiple curve configurations for ease of placement
- Push/pull handle for steering control
- 1 mm band electrodes


Ordering Information

4 F Quadripolar steerable diagnostic catheter (1 unit per box)

| Reorder Number | Description | Tip Electrode (mm) | Electrode Spacing (mm) | Curve Type | Usable Length (cm) |
|----------------|-------------|--------------------|------------------------|------------|--------------------|
| 81540 | 1104-4-25-M | 1 | 2-5-2 | Medium | 110 |
| 81542 | 1104-4-5-M | 1 | 5 | Medium | 110 |
| 81543 | 1104-4-25-L | 1 | 2-5-2 | Large | 110 |

Required Catheter Connecting Cable – Page 122

| Reorder Number | Model Number | Description | Length (m) |
|----------------|--------------|-----------------------------------|------------|
| 85953 | 1904-SA | 4-Pin Diagnostic Connecting Cable | 1.5 |



Product Highlights

- Variable electrode spacing for pacing and recording
- Multiple curve configurations for ease of placement
- Push/pull handle for steering control
- 1 mm band electrodes

Ordering Information

4 F Decapolar steerable diagnostic catheter (1 unit per box)

| Reorder Number | French Size | Electrode Spacing (mm) | Curve Type | Usable Length (cm) |
|----------------|---------------------|------------------------|------------|--------------------|
| 81530 | 1110-4-2-M | 2 | Medium | 110 |
| 81531 | 1110-4-25-M | 2-5-2 | Medium | 110 |
| 81532 | 1110-4-25-L | 2-5-2 | Large | 110 |
| 81534 | 1110-4-5-L | 5 | Large | 110 |
| 81537 | 1110-4-25-M(SC)(60) | 2-5-2 | M(SC) | 60 |

Note: (60) indicates catheter length in cm

Required Catheter Connecting Cable – Page 122

| Reorder Number | Model Number | Description | Length (m) |
|----------------|--------------|------------------------------------|------------|
| 85954 | 1910-SA | 10-Pin Diagnostic Connecting Cable | 1.5 |
| 85930 | 1910-S | 10-Pin Diagnostic Connecting Cable | 1.5 |
| 85942 | 1910-S | 10-Pin Diagnostic Connecting Cable | 2.5 |

