

# Pacel™ Bipolar Pacing Catheter, Straight

**Flow Directed**

## Product Highlights

- Designed to provide quick, reliable temporary pacing in emergency situation with or without the use of fluoroscopy
- Flexible design for delicate placement of catheter
- Rapid inflation and deflation of balloon
- 10 cm depth markers to enhance placement
- Smooth, hemocompatible polyurethane surface
- Clearly marked terminals
- Shrouded 2 mm pin connector leads

## Ordering Information

401761 Contents: Pacing Catheter, Syringe and Pin Adapters (1 unit per box)

401615 Contents: Pacing Catheter and Syringe (1 unit per box)

401762 Kit Contents: Pacing Catheter, Syringe, Pin Adapters and Needle Peel-Away Introducer (1 unit per box)

401616 Kit Contents: Pacing Catheter, Syringe, Needle Peel-Away Introducer and ECG Adapter (1 unit per box)

Shrouded Pin Reorder Number	Exposed Pin Reorder Number*	French Size	Electrode Spacing (cm)	Usable Length (cm)	Recommended Introducer Size (F)
401761	401615*	5	1	110	6
401762 (Kit)	401616 (Kit)	5	1	110	6



\*Allow up to 3 weeks for delivery.

\*Not available in countries requiring CE Mark.

Product referenced is approved for CE Mark.

## Pacel™ Bipolar Pacing Catheter, Right Heart Curve

### Flow Directed

### Product Highlights

- Designed to provide quick, reliable temporary pacing in emergency situation with or without the use of fluoroscopy
- Flexible design for delicate placement of catheter
- Rapid inflation and deflation of balloon
- 10 cm depth markers to enhance placement
- Smooth, hemocompatible polyurethane surface
- Clearly marked terminals
- Shrouded 2 mm pin connector leads

### Ordering Information

401763 Contents: Pacing Catheter, Syringe and Pin Adapters (1 unit per box)

401623 Contents: Pacing Catheter and Syringe (1 unit per box)

401764 Kit Contents: Pacing Catheter, Syringe, Pin Adapters and Needle Peel-Away Introducer (1 unit per box)

401624 Kit Contents: Pacing Catheter, Syringe, Needle Peel-Away Introducer and ECG Adapter (1 unit per box)

Shrouded Pin Reorder Number	Exposed Pin Reorder Number*	French Size	Electrode Spacing (cm)	Usable Length (cm)	Recommended Introducer Size (F)
401763	401623*	5	1	110	6
401764 (Kit)	401624 (Kit)	5	1	110	6



\*Allow up to 3 weeks for delivery.

\*Not available in countries requiring CE Mark.

Product referenced is approved for CE Mark.